

UAE

Email: mathews@generaltechuae.com

Tel: +971 (0) 2 550 7702 Fax: +971 (0) 2 550 7706 SAUDI

Mathews@generaltechsaudi.com

Tel: +966 (0) 3 834 4264 Fax: +966 3 834 4538



P.O Box 25898 Sharjah

United Arab Emirates

Apera Instruments AI1102 | Standard Calibration Buffer Solution | Accuracy ±0.01 pH

Color-coded and calibrated to NIST-certified pH guidelines, the Apera Instruments AI1102 pH 4.00 Calibration Buffer Solution is within 0.01 pH at 25°C (77°F). It is perfect for calibrating various pH meters and testers. A temperature/pH reference chart is included on the bottle label for reference since temperature has an impact on pH level. Every solution has a two-year shelf life, is double-sealed, and is exempt from DOT shipping regulations. The bottle has an expiration date printed on it. MSDS documents are available for download. On request, a certificate of analysis can be obtained.

Specifications Details

 Phone
 : 06-5436933/06-5436733
 E- Mail
 : gentech@eim.ae

 Fax
 : 06-5437077
 Website: www.generaltech.ae